



**ThermoFisher**  
S C I E N T I F I C

## Flexible Panel-Based Microbial Detection using Spatial Multiplexing on Nanofluidic qPCR Platform

Nitin Puri, PhD

Associate Director, Infectious Disease

Oct 29<sup>th</sup>, 2018

The world leader in serving science

# Applied Biosystems - Microbial Detection Solutions


## End to end solution approach for supporting pathogen testing

Organism Name	Target	Assay Name ID
Staphylococcus aureus	StaphSA001	StaphSA001_01
Staphylococcus aureus	StaphSA002	StaphSA002_01
Staphylococcus aureus	StaphSA003	StaphSA003_01
Staphylococcus aureus	StaphSA004	StaphSA004_01
Staphylococcus aureus	StaphSA005	StaphSA005_01
Staphylococcus aureus	StaphSA006	StaphSA006_01
Staphylococcus aureus	StaphSA007	StaphSA007_01
Staphylococcus aureus	StaphSA008	StaphSA008_01
Staphylococcus aureus	StaphSA009	StaphSA009_01
Staphylococcus aureus	StaphSA010	StaphSA010_01
Staphylococcus aureus	StaphSA011	StaphSA011_01
Staphylococcus aureus	StaphSA012	StaphSA012_01
Staphylococcus aureus	StaphSA013	StaphSA013_01
Staphylococcus aureus	StaphSA014	StaphSA014_01
Staphylococcus aureus	StaphSA015	StaphSA015_01
Staphylococcus aureus	StaphSA016	StaphSA016_01
Staphylococcus aureus	StaphSA017	StaphSA017_01
Staphylococcus aureus	StaphSA018	StaphSA018_01
Staphylococcus aureus	StaphSA019	StaphSA019_01
Staphylococcus aureus	StaphSA020	StaphSA020_01

**COVERAGE**  
Allows to select clinical research studies and population demographics relevant genomic targets



**INSTRUMENTATION**  
Combining flexible throughput capabilities with a streamlined workflow



**PLATFORM**  
Form factors to suit different sample number, target number and throughput needs




**CONSUMABLES**  
Supporting workflow reagents including controls, sample prep, master mixes



**SERVICE & SUPPORT**  
Instrument OQ/IQ/PV, training and support your analytical validation



# Applied Biosystems Solution for Real-Time PCR Detection of Microbiota



Taqman<sup>®</sup> **OpenArray<sup>®</sup>** (OA) Platform offers a a pan-molecular detection method for quantifiable and comparable results of all tested microbes, thereby reducing the variability associated with multiple conventional methods.



**High Throughput**  
Up to 800  
Samples a day



**Performance Proficiency**  
High Sensitivity and  
Specificity



**Flexibility & freedom**  
Selecting only  
what you  
want to test for to  
suit your needs




**Extensive Coverage**  
Range of  
commensal and  
pathogenic  
microbes



**Unrivalled support**  
Validation services  
& 24/7 support  
available


# OpenArray™ Anatomy




**Flexibility to customize your own panels  
for number of samples and targets**



# Flexibility for Customization for Targets & Samples




- Spatial Multiplexing allows for 56 Molecular targets per subarray
- Selection of assays for commensal and pathogenic microbes, host factors, antibiotic resistance genes



Broncho alveolar Lavage



Nasopharyngeal Aspirate



Nasopharyngeal Swab

# Blood-borne pathogen study using OA (Grigorenko et al, 2017)

Study from Robert Duncan et. al. at Center for Biologics Evaluation and Research, FDA, simultaneous detection and discrimination of **17 viral pathogens** in human plasma samples and **13 bacterial and protozoan pathogens** in human blood samples on the **OpenArray platform**. The custom plate was tested for specificity and analytical sensitivity with purified nucleic acids from each pathogen and with pathogen-spiked human blood and plasma samples.

The Journal of Molecular Diagnostics, Vol. 19, No. 4, July 2017



## Highly Multiplex Real-Time PCR-Based Screening for Blood-Borne Pathogens on an OpenArray Platform



Elena Grigorenko,\* Carolyn Fisher,<sup>1</sup> Sunali Patel,<sup>1</sup> Valerie Winkelman,<sup>1</sup> Philip Williamson,<sup>1</sup> Caren Chancey,<sup>1</sup> Germán Afnez,<sup>1</sup> María Ríos,<sup>1</sup> Victoria Majam,<sup>1</sup> Sanjai Kumar,<sup>1</sup> and Robert Duncan<sup>1</sup>

From the Department of Research and Development,\* Diagnostics Laboratories, LLC, Huntsville, Alabama; the Division of Emerging and Transfusion Transmitted Diseases,<sup>1</sup> Center for Biologics Evaluation and Research, Food and Drug Administration, Silver Spring, Maryland; the Genetic Sciences Division,<sup>1</sup> Thermo Fisher Scientific, South San Francisco, California; and Research and Development Support<sup>1</sup> and Operations and Scientific Affairs,<sup>1</sup> Creative Testing Solutions, Tempe, Arizona

Accepted for publication March 28, 2017.

Address correspondence to Robert Duncan, Ph.D., 10903 New Hampshire Ave., Silver Spring, MD 20993. E-mail: rduncan@fda.hhs.gov.

Molecular diagnostics are increasingly used in the blood bank industry. A device that can combine simultaneous detection of multiple targets with the flexibility of inclusion of emerging pathogens is desirable for testing blood products. A highly multiplexed blood-borne pathogen panel (BBPP) using dual-label probe chemistry (TaqMan assays) was developed for simultaneous detection and discrimination of 17 viral pathogens in human plasma samples and 13 bacterial and protozoan pathogens in human blood samples on the OpenArray platform. The custom BBPP OpenArray plate was tested for specificity and analytical sensitivity with purified nucleic acids from each pathogen and with pathogen-spiked human blood and plasma samples. The results of analytical validation of known samples yielded decision trees for identification of coded samples; pathogens spiked in human plasma or whole blood. Results from coded samples demonstrated no false positives among the plasma or whole blood specimens. Samples not detected were at the lower limit of the detectable range or qualified for retesting as indeterminate. Further demonstration of the performance of the BBPP OpenArray was achieved with clinical samples from a blood donor testing organization. Ninety-five percent of virus-positive samples were correctly identified. These results show that a high-throughput OpenArray PCR platform can be expanded and adapted for higher discrimination and newly emerging agents, enabling consideration for development as a next-generation device for testing blood products. (*J Mol Diagn* 2017; 19: 549–560; <http://dx.doi.org/10.1016/j.jmoldx.2017.03.004>)

The quality and safety of blood products used in transfusions are of significant public concern. The number of emerging pathogens that affect blood safety has steadily increased in the recent decades and has triggered a rapid development of molecular tests for detection of blood-borne pathogens such as West Nile virus,<sup>1</sup> dengue virus (DENV),<sup>2</sup> and *Babesia microti*<sup>3</sup> to name a few. Some highly virulent pathogens may have a low prevalence rate and/or be restricted seasonally or geographically, suggesting they may not be of concern for blood safety. However, the impact of transfusion-transmitted infection of such agents can have fatal consequences in highly vulnerable populations such as

Supported by the Food and Drug Administration Medical Countermeasures Initiative (R.D.).

Disclosures: E.G. is an employee of Diagnostics Laboratories, S.P. is an employee of Thermo Fisher Scientific, the manufacturer of the OpenArray; V.W. and P.W. are employees of Creative Testing Solutions, the source of blood donor specimens.

The comments in this publication are an informal communication and represent the authors' own best judgment. These comments do not bind or obligate the Food and Drug Administration (FDA). The device described in this publication is for research use only and has not been formally cleared or approved by FDA for the uses discussed herein.

All samples and Thermo Fisher Scientific (Applied Biosystems) reagents and materials are for research use only, not for use in diagnostic procedures. Current address of G.A., Sanofi Pasteur, Swiftwater, PA.

**Table 3** Analytical Sensitivity Study Results of Whole Blood Panel

Pathogen	Cells/mL	Assay 1			Assay 2				
		Name	Mean Cq	95% CI	TH Pos/Total	Name	Mean Cq	95% CI	TH Pos/Total
<i>Babesia microti</i>	1000	BAB-1	22.90	22.67–23.13	18/18				
<i>B. microti</i>	100	BAB-1	26.31	25.8–26.82	18/18				
<i>Leishmania braziliensis</i>	1000	LCHAG	26.31	25.66–26.96	18/18	LEI-1	23.58	±0.29	18/18
<i>L. braziliensis</i>	100	LCHAG	24.08	23.81–24.35	18/18	LEI-1	21.84	±0.24	18/18
<i>L. donovani</i>	1000	LINF	23.21	22.77–23.65	18/18	LTRO	22.18	±0.40	18/18
<i>L. donovani</i>	100	LINF	26.70	26.12–27.28	18/18	LTRO	25.82	±0.49	18/18
<i>L. infantum</i>	1000	LCHAG	21.39	21.15–21.63	18/18	LINF	24.70	±0.36	18/18
<i>L. infantum</i>	100	LCHAG	25.06	24.59–25.53	18/18	LINF	28.68	±0.46	18/18
<i>L. major</i>	1000	LMAJ	26.29	25.71–26.87	18/18	LTRO	28.58	±0.57	18/18
<i>L. major</i>	100	LMAJ	25.78	25.4–26.16	18/18	LTRO	28.62	±0.45	18/18
<i>L. mexicana</i>	1000	LMEX	24.18	24.01–24.35	18/18	LMEX-1	27.49	±0.48	17/18
<i>L. mexicana</i>	100	LMEX	28.25	27.96–28.54	12/18	LEI-1	26.03	±0.41	11/18
<i>L. tropica</i>	1000	LTRO	23.97	23.67–24.27	18/18	LEI-2	24.56	±0.31	18/18
<i>L. tropica</i>	100	LTRO	26.90	26.37–27.43	18/18	LEI-2	28.48	±0.50	17/18
<i>Plasmodium falciparum</i>	1000	PLA-3	21.58	21.19–21.97	17/18	PLA-2	19.95	±0.31	18/18
<i>P. falciparum</i>	100	PLA-3	25.12	24.5–25.74	18/18	PLA-2	23.06	±0.35	18/18
<i>Plasmodium vivax</i>	1000	PLA-4	23.49	23.11–23.87	18/18	PLA-2	18.88	±0.36	18/18
<i>P. vivax</i>	100	PLA-4	25.58	24.99–26.17	18/18	PLA-2	21.37	±0.32	18/18
<i>Trypanosoma cruzi</i>	1000	TCF-1	13.47	12.95–13.99	18/18	TCF-2	13.40	±0.53	18/18
<i>T. cruzi</i>	100	TCF-1	17.35	16.86–17.84	18/18	TCF-2	17.20	±0.46	18/18
<i>Escherichia coli</i>	1000	GNEG-1	24.39	24.1–24.68	17/18	GNEG-2	24.89	±0.36	16/18
<i>E. coli</i>	100	GNEG-1	27.31	26.81–27.81	18/18	GNEG-2	27.40	±0.36	18/18
<i>Yersinia enterocolitica</i>	1000	GNEG-1	25.21	24.87–25.55	18/18	GNEG-2	25.71	±0.27	18/18
<i>Y. enterocolitica</i>	100	GNEG-1	28.98	28.43–29.53	15/18	GNEG-2	29.24	±0.56	16/18
<i>Staphylococcus aureus</i>	10,000	STAU	25.11	24.73–25.49	18/18				
<i>S. aureus</i>	1000	STAU	27.67	26.91–28.43	5/18				
Neg control blood		LEI-1	28.66	28.26–29.06	15/27	GNEG-1	29.53	±0.37	4/27
No-template control		GNEG-1	29.38	29.03–29.73	10/36	GNEG-2	29.64	±0.45	9/36

# Sample-to-Answer Workflow for UTM

Up to 3 runs per day

Total Turnaround Time  
5 Hours / Run



Service and support available throughout the entire process

Sample type: urine



Sample prep

<2 Hrs



Thermo Scientific™ KingFisher™ Flex system with Applied Biosystems™ MagMAX kits

Sample loading

20 min



Applied Biosystems™ OpenArray™ AccuFill™ System for sample loading


Reporting



Data Analysis  
In-house LIS or CLS




Analysis



Run real-time PCR

2 hr




4 plates x 12 per day

UTI Pathogen Detection Assays

# Microbiome Applications Areas in use on OpenArray


**Urogenital**



**Gastrointestinal**



**Respiratory**



**Antibiotic Resistance**




**Wound**



**Foot Fungus**



**Tuberculosis**




**CNS**







**QuantStudio**



**OpenArray**

- OpenArray® and QuantStudio offer a low cost, customized and proficient microbial detection solution system
- Applied Biosystems offers a tool box with supporting reagents including sample prep and customized controls'

Reach out to [nitin.puri@thermofisher.com](mailto:nitin.puri@thermofisher.com) for potential interest / collaborations



Thank you

**ThermoFisher**  
S C I E N T I F I C

**For Research Use Only. Not for use in diagnostic procedures.** © 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. **COL15382 0218**

The world leader in serving science